

ADAPTATION OF COMPRESSION CLASSES

Following a company review phase, it was decided to reclassify the compressions detected on the limb according to the cataloging described in the German RAL-GZ 387 standard, currently considered the most authoritative.

The compressive force exerted on the products has remained unchanged, the detection system has been updated and takes place through new certified equipment.

Consequently, the compression values reported by the new system differ slightly from the data previously indicated although the products remain perfectly the same in terms of compression.

The product packages will no longer show the mm Hg indication on the label but only the Compression Class (Ccl 1, Ccl 2, Ccl 3).

The values detected by the new certified equipment pass:

- for Ccl 1 from 20-30 mm Hg to 18-21 mm Hg**
- for Ccl 2 from 30-40 mm Hg to 23-32 mm Hg**
- for Ccl 3 from 40-50 mm Hg to 34-46 mm Hg**